


Wi-Fi Simultaneous Image on Secondary Display

- › Aitraq Max View 10.4" Display
- › Phone, Tablet and PC

Aitraq
MaxView



Phone

Tablet

PC

Additional wire connection to Aitraq MaxView Display and to most HDMI monitors

Aitraq
+Plus

One VL fits **All** blade types & sizes



MacIntosh



Hyper Angulated



Channeled

Blades made from
RECYCLED Polycarbonate
Your **ECO** Choice



www.airtraq.com

USA
Airtraq LLC
Phone: 1-877-624-7929
E-mail: info.usa@airtraq.com

EU, Asia & RoW
Prodol Meditec S.A.
Phone: +34 944 804 690
E-mail: info@airtraq.com

Aitraq
+Plus

All In One Eco
Video laryngoscope



VERSATILE
AFFORDABLE

SUSTAINABLE
SAFE

Intubation from Any Position in Any Scenario



For Routine & Difficult Airways



Perfectly balanced



Full hermetic blades
with integrated handle



3.5" touch screen.
User friendly navigation



Video Laryngoscope into
Charging Dock

Specifications

Reference	A600
Description	Video Laryngoscope
Screen	3.5" IPS. Resolution 640 x 480
Touch Panel	Viewing angle = 85°
Video Sensor	CMOS. Resolution: 800 x 800
Battery Capacity	180 min. Rechargeable
Wi-Fi	Range: 5 m (approx.)
Internal Memory	16 Gb. 10 hours (approx.)
IP Protection	IP6x (IEC 60529)
Weight	195 grams
Auto-Recording	Yes. Selected by user
Playback	Yes. Full functionality
Hybrid Intubation	Yes. Manual control of VL light

Single Use **ECO** Blades



REF.	DESCRIPTION	SIZE	
A6M4	MacIntosh # 4	4-Adult	> All 3 blade geometries > From neonates to adults > Made from ISCC Certified Recycled Polycarbonate
A6M3	MacIntosh # 3	3-Small Adult	
A6M2	MacIntosh # 2	2-Pediatric	
A6M1	MacIntosh # 1	1-Neonate	
A6HA	Hyper Angulated "H"	Adult	> One-piece handle and blade minimizes cross-contamination > Low profile > Hermetic
A6C3	Channeled Regular	3-Adult	
A6C2	Channeled Small	2-Small Adult	
A6C1	Channeled Pediatric	1-Pediatric	
A6C0	Channeled Infant	0-Neonate	