

A390 Wi-Fi Camera



iOS / Android
MAC / Windows

- > Video & Snapshot recording
- > Auto recording
- > Full on-device playback
- > Easy video transfer to PC



USA
Airtraq LLC
Phone: 1-877-624-7929
E-mail: info.usa@airtraq.com

www.airtraq.com

EU, Asia & RoW
Prodol Meditec S.A.
Phone: +34 944 804 690
E-mail: info@airtraq.com

Airtraq sp

Channeled Video laryngoscope



**Fully Disposable
One Piece
Ready to Use**

Simple ♦ Effective ♦ Affordable

Airtraq *sp*

Available in all sizes
from Neonate to Adult

■ Infant	■ Pediatric	■ Nasal	■ Double lumen
Ref: A041	Ref: A031	Ref: A061	Ref: A071
2.5-3.5 ETT	4.0-5.5 ETT	-	28-41 Fr
M.O.*: 11 mm	M.O.*: 12 mm	M.O.*: 15 mm	M.O.*: 18 mm

(*) M.O.: Minimum Mouth Opening

■ Adult regular	■ Adult small
Ref: A011	Ref: A021
7.0-8.5 ETT	6.0-7.5 ETT
M.O.*: 16 mm	M.O.*: 15 mm



Specially effective for:

- › Cervical immobilization
- › Obese
- › Awake
- › In combination with a fiberscope
- › Frontal or lateral intubations
- › Poor dentition

Eyecup connects to
endoscopy cameras



How Channeled VL simplifies intubation

- › Airtraq's channel simplifies ETT insertion reducing intubation time.
- › No stylet needed. Avoids potential injuries. Eliminates stylet costs.
- › 90° shape minimizes hyperextension and reduces force required.
- › Reduces hemodynamic reactivity.
- › Facilitates intubation from any position (face to face, etc).
- › High intubation success rate in difficult airways.
- › Easier to learn than other video laryngoscopes.

90° shape
minimizes
force
exertion

Channel
simplifies
ETT
insertion

Available
when and
where you
need it