

How Airtraq simplifies Video Intubation

- Airtraq's channel simplifies ETT insertion reducing intubation time.
- No stylet needed. Avoids potential injuries. Eliminates stylet costs.
- 90° shape minimizes hyperextension and reduces force required.
- Reduces hemodynamic reactivity.
- Facilitates intubation from any position (face to face, etc.).
- High intubation success rate in difficult airways.
- Easier to learn than other video laryngoscopes.

Multiple viewing options

90° channel simplifies intubation

Minimum capital investment allows availability in every location

Visualization Options

for any Airtraq model



Airtraq LLC
2414 Lawton Ln.
Rowlet, TX 75089, USA
Phone: 1-877-624-7929
E-mail: info.usa@airtraq.com

Airtraq
GUIDED VIDEO INTUBATION
www.airtraq.com

Prodol Meditec S.A
Muelle T. Olavarrí, 5
48930 Las Arenas, Spain
Phone: +34 944 804 690
E-mail: info@airtraq.com

Airtraq
GUIDED VIDEO INTUBATION



Simple, effective,
affordable video
laryngoscopy

Two complementary models for one complete solution



AIRTRAQ AVANT


Reusable Optics + Disposable Blades

Minimum cost per intubation.
VL for frequent use.

AIRTRAQ SP

Fully Disposable

One piece. Ready to use.
VL for remote locations and emergencies.



Airtraq improves routine intubations and it is specially effective for:

- Cervical immobilization
- Obese
- Awake
- In combination with a fiberscope
- Frontal or lateral intubations
- Poor dentition

Airtraq sp

Adults
Pediatrics
Double Lumen
Nasal (no channel)



Adult regular
7.0-8.5 ETT
M.O.*: 16 mm

Adult small
6.0-7.5 ETT
M.O.*: 15 mm

Double lumen
28-41 Fr
M.O.*: 18 mm

Nasal
M.O.*: 15 mm



Pediatric
4.0-5.5 ETT
M.O.*: 12 mm

Infant
2.5-3.5 ETT
M.O.*: 11 mm

NEW

Wi-Fi Camera

Airtraq CAM

iOS / Android
Mac / Windows



Phone adapter



iOS / Android



Adult regular
7.0-8.5 ETT
M.O.*: 17 mm

Adult small
6.0-7.5 ETT
M.O.*: 17 mm



Reusable Optics

Disposable blades

Docking

Airtraq Avant

* Mouth opening