

New Camera A390 vs. Previous A360

- ✓ **Weight**: 110 grams. **50% less** than A360.
- ✓ **Size**: **40% less** volume than A360.
- ✓ **Screen size**: remains the same generous 2.8 inches.
- ✓ **Impact Resistance**: **Integrated bumpers** and solid one piece construction without moving parts.
- ✓ **Attachment to Airtraq**: Secure friction lock. No additional latch.
- ✓ **Boot up time**: **4 secs.** vs 45 secs.
- ✓ **Battery time without Wi-Fi**: **240 min** vs. 120 min.
- ✓ **Battery time with Wi-Fi**: 180 min vs. 120 min.
- ✓ **Image for Pediatric mode**: improved focus
- ✓ **User Settings**: App on PC for easier configuration vs. Camera screen.
- ✓ **Software updates**: App on PC vs. complex manual procedure.
- ✓ **Docking**: Easier connection for charging.
- ✓ **Silicon Case**: 1 piece custom designed for A390 provides better protection for abusive environments.
- ✓ **Semi-Rigid ballistic nylon Storage Case**: Smaller, lighter, more durable.
- ✓ **Neoprene camera bag**: Additional transport protection.
- ✓ **Reduced cost to facilitate availability in every location.**